



TWYDIL® LEG GEL

Horse care product
Gel for tendons and ligaments of horses.

Composition

Iodine	1.0 g	Isopropyl Alcohol	7.5 g
Potassium iodide	0.7 g	Polyethylene Glycol 300	90.8 g

Instructions for use

Apply a thick layer twice a day to the horse's legs without rubbing. Cover with a bandage where possible. Blood circulation will be stimulated locally, leading to a natural regeneration.

For external use only.
Keep out of reach of children.
Batch number and expiry date are indicated on the jar.

Packaging

Jar of 250 g net
Carton of 12 jars of 250 g net

Producer

Pavesco AG, Elisabethenstr. 54, CH-4010 Basel
www.twydil.com

(ref. 05.2015)